

MedCAD Digital Library

Quality Assurance Checklist

To be completed by the Regional AM Hub prior to dispatch of any fabricated component

Document Reference: MedCAD-QA-004 | Version: 1.0 | Classification: DRAFT

Order Reference	
Component Reference	
Component Description	
Requesting NHS Trust	
Regional AM Hub	
Fabrication Date	
Operator Name & Signature	
CAD Datapack Version	
Machine / Printer ID	
Material Batch Reference	

IMPORTANT: All checklist items must be completed before this component is cleared for dispatch. Any item recorded as 'Fail' must be investigated and resolved. No component may be dispatched until all applicable items are recorded as 'Pass' or 'N/A', and this form has been signed off by the QA Authoriser.

Quality Assurance Checklist

Check Item	Acceptance Criteria / Notes	Reference	Result
Section 1 — File & Data Verification			
CAD Datapack version confirmed	Fabrication file version matches the current approved version in the MedCAD library.	QA-1.1	Pass / Fail / N/A
File integrity check passed	Hash or checksum of downloaded file matches Platform-issued value. No corruption detected.	QA-1.2	Pass / Fail / N/A
OEM authorisation status verified	Component is listed as 'OEM Authorised' and not under any hold or suspension in the Platform.	QA-1.3	Pass / Fail / N/A
Permitted scope of use reviewed	Operator has read and understood any OEM-specified restrictions on application or duration.	QA-1.4	Pass / Fail / N/A
Section 2 — Material Verification			
Material grade confirmed	Material used matches OEM-	QA-2.1	Pass /

	specified grade exactly. No substitutions made.		Fail / N/A
Material batch certificate reviewed	Certificate of conformity for material batch is on file and within validity period.	QA-2.2	Pass / Fail / N/A
Material storage conditions verified	Material stored in accordance with supplier requirements prior to use.	QA-2.3	Pass / Fail / N/A
Section 3 — Fabrication Process			
Process parameters verified	Layer height, speed, temperature, infill, and orientation match OEM-approved build parameters.	QA-3.1	Pass / Fail / N/A
Machine calibration current	Printer/machine has been calibrated within the required maintenance interval.	QA-3.2	Pass / Fail / N/A
Build log retained	Complete machine build log exported and stored against this order reference.	QA-3.3	Pass / Fail / N/A
Post-processing completed	All required post-processing steps (support removal, surface finishing, heat treatment, etc.) completed per specification.	QA-3.4	Pass / Fail / N/A
Section 4 — Dimensional & Visual Inspection			
Critical dimensions verified	All critical dimensions measured and within OEM-specified tolerances. Measurement record attached.	QA-4.1	Pass / Fail / N/A
Surface finish acceptable	Surface finish meets OEM specification. No visible defects, delamination, or porosity.	QA-4.2	Pass / Fail / N/A
Thread and mating features checked	All threads, bores, and mating surfaces conform to specification.	QA-4.3	Pass / Fail / N/A
Component free of contamination	Part clean and free from debris, oil, or process residue.	QA-4.4	Pass / Fail / N/A
Section 5 — Traceability and Documentation			
Unique part identifier applied	Permanent or tamper-evident identifier applied to component linking it to this QA record.	QA-5.1	Pass / Fail / N/A
QA record complete	All fields in this checklist completed and signed. All supporting records attached.	QA-5.2	Pass / Fail / N/A
Digital Thread updated	Fabrication event, QA results, and operator details logged in Platform audit trail.	QA-5.3	Pass / Fail / N/A
Component marked as temporary replacement	Packaging or accompanying documentation clearly states component is a temporary replacement pending OEM supply.	QA-5.4	Pass / Fail / N/A

Overall Disposition

Disposition Decision	☐ CLEARED FOR DISPATCH — REQUIRES REVIEW	☐ HOLD ☐ REJECTED
Reason (if Hold or Rejected)		

QA Authoriser Sign-Off

QA Authoriser Name	Signature	Date	Hub Reference

This QA Checklist and all supporting documentation must be retained for a minimum of seven (7) years from the date of fabrication. Records must be stored in the MedCAD Platform audit trail and, where applicable, in the Regional AM Hub's quality management system.